

Biochemistry Incubator

Features:

- * BJPX series: LED display.
- * BJPX-II series: LCD display
- * Polished stainless steel chamber.
- * Round angle structure, easy to clean.
- * Universal castors for easy movement.
- * CFC-free eco-friendly refrigerant.
- * Microprocessor PID temperature control, accurate and reliable.
- * Power socket integrated inside for additional small equipment.



BJPX-B

BJPX-BII

BJPX-BI

Safety Function:

- * Temperature protection threshold value can be set.
- * Auto-off protection and alarm for abnormal situation.
- * Parameters can be saved automatically when power off or system error.



USB Adaptor



RS485 Interface

Optional Accessories:

- * Printer (for type II series)
- * USB adaptor
- * Independent temperature protector
- * Test port

Technical Parameters:

Model	BJPX-B80	BJPX-B150	BJPX-B200	BJPX-B250I	BJPX-B250	BJPX-B300	BJPX-B 400
	BJPX-B80II	BJPX-B150II	BJPX-B200II	BJPX-B250I	BJPX-B250II	BJPX-B300II	BJPX-B400II
Capacity	80L	150L	200L	250L	250L	300L	400L
Temp. Range	5~50°C(type II 0~65°C); B250I: 0-60°C						
Temp. Fluctuation	±1°C (type II: heating: +/-0.5°C); B250I: ±0.5°C						
Temp. Uniformity	±2°C(type II: ±1); B250I: ±1°C						
Time Range	0~9999min						
Refrigerant	R134a (CFC free, eco-friendly)						
Standard Accessory	BOD socket						
Consumption	180W	250W	300W	300W	300W	400W	550W
Power Supply	AC220V±10%, 50Hz						
Internal Size(mm)	400*370*550	450*420*850	450*450*1000	510*430*1085	480*490*1070	520*500*1170	580*540*1270
External Size(mm)	540*570*1070	590*620*1370	590*640*1520	595*604*1635	620*680*1590	660*690*1690	720*740*1790
Package Size(mm)	700*730*1260	750*780*1560	750*800*1710	680*690*1785	780*840*1780	820*850*1880	880*900*1980
Gross Weight(kg)	85	120	138	120	147	157	170

Biochemistry Incubator

Features:

- * LCD display Controller.
- * Illumination in the chamber.
- * Polished stainless steel chamber.
- * Round angle structure, easy to clean.
- * Universal Castors for easy movement.
- * CFC-free environmental friendly refrigerant.
- * Microprocessor PID temperature control: accurate and reliable.
- * Power socket integrated inside for additional small equipment.



BJPX-BIII



Micro Printer

Safety Function:

- * Temperature protection threshold value can be set.
- * Auto-off protection and alarm for abnormal situation.
- * Parameters can be saved automatically when power off or system error.

Optional Accessories:

- * Independent temperature protector
- * Test port

Technical Parameters:

Model	BJPX-B150III	BJPX-B200III	BJPX-B250III	BJPX-B300III	BJPX-B400III
Capacity	150L	200L	250L	300L	400L
Temp. Range	0~65°C				
Temp. Fluctuation	±1°C				
Temp. Uniformity	±2°C				
Shelves	2 pcs/adjustable				
Time Range	0~99h 59min				
Refrigerant	R134a (CFC free, eco-friendly)				
Working Type	Continuous				
Standard Accessory	BOD socket; micro printer				
Consumption	300W	350W	450W	500W	600W
Power Supply	AC110/220V±10%, 50/60Hz				
Internal Size(mm)	590*590*1420	590*640*1560	620*680*1630	660*690*1730	720*740*1830
External Size(mm)	450*420*850	450*450*1000	480*490*1070	520*500*1170	580*540*1270
Packing Size(mm)	730*730*1600	750*800*1750	820*770*1820	840*800*1920	870*870*2020
Gross Weight(kg)	120	138	147	157	170